

INSTRUCTIONS FOR USE

I.V.CANNULA WITH SMALL WINGS AND WITHOUT INJECTION PORT

MATERIAL USED :

-PTFE/PUR/FEP,PP,HDPE,POM , Silicon Rubber & Stainless Steel

INDICATIONS :

- Infusion of I.V. Solutions
- To maintain hydration and/or correct dehydration if patient is unable to take sufficient volume of oral fluid.

CONTRAINDICATIONS:

- Product should not to be used in patients with known hypersensitivity to any of the materials used.
- Administration of high viscous fluids.
- Large blood transfusion.

INSTRUCTION FOR USE :

- Carefully select and aseptically prepare the site.
- Select suitable size of I.V. Cannula & inspect visually to ascertain that package is intact.
- Remove cannula from sterile packing and remove its needle cover.
- Perform veinpuncture & check for flash back of blood in flash back chamber.
- Advance the catheter into vein and simultaneously withdraw the needle.
- Catheter & Needle is coated with Polydimethylsiloxane.
- The product must not be used in CT Scan and MRI.
- Re-use of this device may change its mechanical or biological features and may cause device failure, allergic reactions of infections.
- This product if contains phthalate should not be used on children, pregnant and nursing woman and other patient at risk.
- Contains universal 6% taper that can be connected to any standard medical device.

NEVER TRY TO REINSERT THE PARTIALLY OR COMPLETELY WITHDRAWN NEEDLE.

- Withdraw the needle completely while pressing the vein just after the tip of catheter into the vein & discard the needle in an appropriate container.
- Connect to the I.V. infusion set line.
- Cover the puncture site with sterile dressing.
- Perform routine monitoring & veinpuncture site maintenance according to medical norms.

WARNING :

- The use of this product is restricted to a qualified doctor or a Paramedic
- Read Instruction before use.

MANUFACTURER DISCLAIMS ANY RESPONSIBILITY FOR POSSIBLE CONSEQUENCES RESULTING FROM IMPROPER USE.

- The product should not be reprocessed.
- Visually inspect and carefully check the product and packaging before use. Improper transport and handling may cause structural and/or functional damage to device or packaging.
- The product is guaranteed sterile & Non-pyrogenic if the package has not been opened or damaged.
- Do not clean or resterilise. For single use only. Discard after use.
- The product should be used immediately after opening the packaging.
- The device must not be re-used after its usage, it can be a potential biohazard.

STORAGE CONDITION :

- Do not expose to heat or direct sunlight.

Disposafe Health and Life Care Limited
Plot No. 1 & 2, Phase-II, Sector-59, Ballabhgarh,
Faridabad-121004, Haryana (INDIA)
E-mail : info@disposafe.com, Web : www.disposafe.com
Customer Care No. 91-129-4170201

EC REP

CEpartner4U B.V.,
ESDOORNLAAN 13, 3951 DB MAARN
THE NETHERLANDS
E-mail: office@cepartner4u.com

REF

Reference No.

LOT

Lot No./Batch No.



Manufacturer



Date of Manufacturing



Use By / Expiry Date



Sterilised By Ethylene Oxide
Single Sterile Barrier System

MD

Medical Device



Do Not Reuse



Do Not Resterilize



Non-Pyrogenic



Do Not Use If Package is Damage



Caution



Consult Instructions For Use

EC REP

Authorised Representative
In The European Community



Temperature Limit


10°C

45°C


① Catheter Insertion




② Catheter is pushed inside the vein & needle is withdrawn



③ Sometimes due to no flash back paramedic re-insert the needle



④ Due to re-insertion the catheter can cut & moves in vein



CE 0123