

DISPOWING

I.V. Cannula with Wings and Without Injection Port

MATERIALS USED:

- PUR/FEP, PP, HDPE, POM & Stainless Steel.

COATING MATERIAL:

- Polydimethylsiloxane

INDICATIONS:

- Blood transfusion or Infusion of I.V. solutions suitable for administration via peripheral veins.
- Intermittent intravenous Drug administration.
- Prophylactic creation of a secure venous access in patients who may require urgent intravenous drug administration, in particular prior to diagnostic or therapeutic procedures.
- This device may also be intended for blood sampling.
- This product is suitable, to be used within CT scan and MRI procedures.

CONTRAINDICATIONS :

- Product should not be used in patient with known hypersensitivity to any of the materials used, including coating material.
- Administration of highly viscous fluids.
- Large volume blood transfusion.

INSTRUCTIONS FOR USE:

- Carefully select and aseptically prepare the insertion site.
- Select suitable size of I.V. Cannula & inspect visually to ascertain that package is intact.
- Remove cannula from sterile packing and twist the needle cover to remove it.
- Grip the cannula from projection provided on catheter body and hub.
- Perform venipuncture & check for flash back of blood in flash-back chamber.
- Advance the catheter into vein, while withdrawing the needle. **NEVER TRY TO REINSERT A PARTIALLY OR COMPLETELY WITHDRAWN NEEDLE INTO THE CATHETER, AS THE LATTER MAY BE CUT OFF, LEADING TO CATHETER EMBOLISM.**
- Withdraw the needle completely while pressing the vein just after tip of catheter to prevent spillage of blood & discard the used needle in a appropriate container.
- Connect to infusion line and cover puncture site with a sterile dressing.
- Perform routine monitoring of the insertion as per the country or hospital protocol.

DURATION OF USE:

- Change according to CDC Guidelines and / or Hospital or

Disposafe Health and Life Care Limited

Plot No. 1 & 2, Phase-II, Sector-59, Ballabhgarh,
Faridabad-121004, Haryana, India
E-mail : info@disposafe.com, Web : www.disposafe.com
Customer Care No. +91-129-4170201

INSTRUCTIONS FOR USE

Institutional protocols.

- Use specialized infusion teams for insertion and monitoring for better patient outcomes.
- The device should be removed in the event of local/systemic symptoms of infection.

WARNINGS:

- The use of this product is restricted to a qualified doctor or a Paramedic.
- The product should be used according to the instructions for use, read the instructions carefully before use.
- The product should not be used for experiment.
- Used product may have potential to biohazards. Handle and dispose of in accordance with accepted medical practice and applicable local, state and country laws and regulations.
- Do not expose to heat or direct sunlight.

DISPOSAFE DISCLAIMS ANY RESPONSIBILITY FOR POSSIBLE CONSEQUENCES RESULTING FROM IMPROPER USE.

- The product should not be re-processed.
- Do not clean or re-sterilize.
- The product and its packaging must be visually inspect before use. Improper transport and handling may cause structural and/or functional damage to device or packaging.
- The product is guaranteed sterile & non-pyrogenic, if the package has not been opened or damaged.
- The product should be used immediately after opening the packaging.
- For single use only, re-use of this device may change its mechanical or biological features and may cause device failure, allergic reactions or infections.
- If there is any change, in expected performance of the device or in case of any malfunction the device should be immediately removed & sent back to supplier/manufacturer for analysis.

CAUTION:

- **DO NOT use scissors OR any sharp tools at or near insertion site.**
- Contamination and/or limited functionality of the device may lead to injury, illness or death of the patient.

EC REP CEpartner4U B.V., ESDOORNLAAN 13,
3951 DB MAARN, THE NETHERLANDS.
E-mail: office@cepartner4u.com

REF

Reference No.

LOT

Lot No./Batch No.



Manufacturer



Date Of Manufacturing



Use By / Expiry Date



Sterilised By Ethylene Oxide
Single Sterile Barrier System



Do Not Reuse



Do Not Resterilize



Non-Pyrogenic



Do Not Use If Package Is Damaged



Caution



Consult Instructions For Use



Authorised Representative
In The European Community



Temperature Limit

MD

Medical Device

