

dispocombi™

ब्लडड्रप स्नमत् स्वबा

INSTRUCTIONS FOR USE

MATERIAL USED :

- H.D.P.E (High Density Poly-Ethylene)

INDICATIONS :

- To be used as a protective stopper

CONTRAINDICATIONS:

- Product should not be used in patients with known hypersensitivity to the material used.

INSTRUCTIONS FOR USE:

- Check visually whether the package is intact.
- Open packaging and remove Injection Stopper from sterile package.
- Connect as desired to the female port of the infusion device in use

WARNING:

- To be use this product by a qualified doctor or a paramedic.
- Read instructions before use.
- The product should be used according to the instruction for use.

Manufacturer disclaims any responsibility for possible consequences resulting from improper use.

- The product should not be reprocessed.
- Visually inspect and carefully check the product and packaging before use.
- Improper transport and handling may cause structural and /or functional damage to device or packaging.
- The product is guaranteed non-toxic, sterile and Pyrogen free if the package is not opened or damaged.
- Do not clean or re-sterilize. For single use only. Discard after use.
- Re-use of this device may change its mechanical or biological features and may cause failure, allergic reactions of infections.



Disposafe Health and Life Care Limited

Plot No. 1 & 2, Phase-II, Sector-59, Ballabhgarh,
Faridabad-121004, Haryana, India

E-mail: info@disposafe.com

Web : www.disposafe.com

Customer Care No. 91-129-4170201



Reference No.



Lot No./Batch No.



Manufacturer



Sterilised By Ethylene Oxide



Date of Manufacturing



Use By / Expiry Date



Do Not Reuse



Do Not Resterilize



Non-Pyrogenic



Do Not Use if
Package is Damaged



Caution



Temp. Limitation



Authorised Representative
in the European Community



Keep away from sunlight



Consult Instructions For Use



CEpartner4U BV, ESDOORNLAAN 13,
3951 DB MAARN, THE NETHERLANDS.
E-mail: office@cepartner4u.com

AW/IFU_1823, Rev:00, Dt: 13-07-2022